

CTS-11141

Acquity® SM Performance Maintenance Kit

*Injector Valve Cartridge (700005236) not included

*Rotor Seal & Pod Repair Kit included

Equivalent to: Waters® 201000174



Technical Data Sheet

1.0 Description

This kit contains all the parts necessary to keep Waters Acquity and Acquity UPLC® Sample Manager systems working at peak performance. It is equivalent to Waters part number 201000174.

2.0 Specifications

2.1 Components

Qty	Part Description
2	2.5mL Wash Syringe (11-3595)
1	Rotor Seal & Pod Repair Kit, Vespel® 3-Groove for Acquity Everflow Injector Pod/Cartridge (700002765) (11-4134)
1	Needle Seal O-Ring, 002 Kalrez® (CTS-11149)
1	Filter, Air, Chassis (CTS-11150)
1	Tube Assy, Inject Outlet (CTS-A11142)
1	Tube Assy, System Outlet, SDS (CTS-A11145)
1	Sample Loop, 10µL (CTS-A11147)
1	Needle Assembly, PEEK™ (CTS-A11151)
1	100µL Sample Metering Syringe, HP (CTS-A11167)

2.2 Materials Specifications

303 SS

304 SS

316 SS

Borosilicate Glass

Glass

Kalrez



Neoprene
PEEK
Polyurethane Foam
PTFE
UHMWPE
Vespel

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.